

CASY | FROM A TO Y: A LIST OF MEASURED CELLS



Is CASY capable of measuring the cells you are interested in?

CASY measures and counts particles in the size from 0.7 µm up to 120 µm.

Three capillaries are used to achieve this huge size range:

- 45 µm capillary for a measuring range of 0.7 – 30 µm (bacteria, small fungi, debris particles...)
- 60 µm capillary for a measuring range of 1.2 – 40 µm (yeast, algae, plankton, sperms, polls...)
- 150 µm capillary for a measuring range of 3.2 – 120 µm (most types of mammalian cells, most types of cell lines, all types of primary cells, big plankton, protozoa...)

A Large number of particles have already been measured with CASY

In order to give you some guidance, this catalogue lists cell types and particles already successfully measured with any CASY – no matter whether it's a TT or a TTC. You might find the organism you are interested in by browsing through the data – but, as this list by far from being complete, you might also find only some quite similar specimen. However, a good starting point for you to set up your specific test samples for CASY measurement.

I have found the organism I'm interested in – what now?

The list on the following pages contains existing setup files from CASY TT and TTC measured by our customers. For comparison purposes, you can request the setup file for a specific sample from us, but please understand that you must not use this file for measurement of your own samples! The listed files are meant to give you an impression of the settings. However, for creating your specific setup, perform a life/dead measurement using [CASYblue](#), as it is described in the [video manual](#) on cellcounting.de or in the [Quick-Manual](#) for CASY TT and CASY TTC. To request a setup file, please send an email to casy@ols-bio.de indicating the serial-number of your instrument and the setup file you request.

This database should grow – you are invited to participate

If you are measuring and counting specimen not listed here, we are happy to receive a setup file! Please send it to casy@ols-bio.de to be included in a growing catalogue!

CASY TT | SPECIMEN MEASURED – SETUP FILES



TT - algae / protozoa

Dinobryon - sweet water; dying in 0.9%NaCl
Mallomonas - sweet water
Ochromonas - sweet water
Paramecium
Selenastrum capricornum

TT - bacteria

Synechococcus PCC7942-10
Synechococcus PCC7942-3

TT - blood cells

blood - human - erythrocytes
blood - human - granulocytes
blood - human - macrophage
blood - human - T-lymphocytes
Blut_Lyse
Platelets

TT - cell lines

25-9 - mouse-mouse - hybridoma; lymphoblast
3T3-NIH - mouse
3T3-NIH - mouse, swiss albino - embryo - fibroblast
A-9 -mouse-connective tissue
A2780
A549 - human - lung; carcinoma-1
A549 - human - lung; carcinoma-2
AG8
ARPE-19 - human - eye, retina; pigmented epithelium
BT-474 - human - mammary gland; ductal carcinoma
C3H-clone8-mouse-embryo fibroblast

Caco-2 - human - colon; adenocarcinoma-1
Caco-2 - human - colon; adenocarcinoma-2
CEpan3b - human - pancreas; stem cells
CHO - clone2
CHO - clone3
CHO - hamster, chinese - ovary cells; transfected
CHO differend clones
CHO transfected with IgG
CHO-K1
COS-7 - monkey, african green - kidney fibroblast
CX-2 - human - colon; adenocarcinoma
Daoy neuroblastoma - human
DU145 - human - brain; carcinoma
HaCaT humane Keratinozyten
HCT-116 - human - colon; colorectal carcinoma
HEK 293 - human - embryonic kidney
HEK 293 - human - kidney; embryonic
HEK 293 - human - kidney; fetal
HEK 293F
HEK293T - human - kidney; highly transfective
HeLa - human - cervix; adenocarcinoma
HeLa - human - cervix; adenocarcinoma
HeLa - human - cervix; adenocarcinoma
HeLa-S3 - human - cervix; adenocarcinoma
Hep-3B- human - liver; hepatocellular carcinoma
HepG2 - human - liver; hepatocellular carcinoma-1
HepG2 - human - liver; hepatocellular carcinoma-2
HKB11
HL-60 - human - promyelocytic leukemia
HL60
HT-29 - human - colon; adenocarcinoma
HUVEC - human - vein; vascular endothelium

CASY TT | SPECIMEN MEASURED – SETUP FILES

IMR-90 - human - lung; fibroblast
 Ishikawa - human - endometrium; carcinoma
 JURKAT - human - T cell; leukemia
 JURKAT - human - T lymphocyte; acute T cell leukemia
 JURKAT - human - T lymphocyte; transgenic
 K-562 - human - pleural effusion; chronic myelogenous leukemia
 K562
 KC - mouse-mouse - hybridoma; lymphoblast
 L5178Y - mouse - lymphoma
 L5178Y - mouse - lymphoma
 L5178Y TK- mouse - lymphoma (clone 3.7.2C)
 LAZ - human - B lymphocyte
 LCL 8664 - monkey, Rhesus - B-lymphocyte; lymphoma
 LN-229 - human - brain; glioblastoma
 LNCaP - human - lymph node; prostate carcinoma
 M1 - human - fibroblast
 McA-3D - rat - liver; hepatoma
 MCF7 - human - adenocarcinoma
 MCF7 - human-adenocarcinoma
 MCF7
 MDBK - cow kidney, cell line
 MLE-12 - mouse, transgenic - lung; epithelial
 PER c6
 PLC-BRF-5 - human - liver; Alexander cells; hepatoma
 Raji - human - B lymphocyte; Burkitts lymphoma
 RAW264-7 - mouse - monocyte; macrophage; AMLV-induced tumor
 RKO-human-colon; carcinoma
 S2 - insect - drosophila; embryo cells
 SCLC-24H - human - lung; small cell carcinoma
 SH-SY5Y - human - bone marrow; (metastasis)
 SH-SY5Y - human - brain; neuroblastoma
 SK-HEP-1 - human - liver; adenocarcinoma
 SL-29 - chicken - embryo; fibroblast
 SW480 - human - colon; adenocarcinoma
 THP-1 - human - acute monocytic leukemia
 TK6 humane Lymphoblasten-2
 TK6 humane Lymphoblasten
 U-266 - human - myeloma
 U-2OS - human - bone; osteosarcoma
 U-373 - human - astrocytes

U-937 - human - macrophage; histiocytic lymphoma
 U937- (Budapest)

TT - parasites

Trypanosoma brucei

TT - primary cells

AF - pig, primary cells
 CCE mouse stem cells-- MEF
 CCE-mouse stem cell and MEF feeder cells
 Embryonic stem cells - human, highly aggregated
 Hepatocytes - rat - after perfusion
 MSC - human - mesenchymal stromal cells, placenta, 2D
 MSC - human - mesenchymal stromal cells, placenta, 3D
 prim - human - chondrocytes
 prim - human - dendritic cells
 prim - human - fibroblast-skin
 prim - human - granulosa -
 prim - human - keratinocytes
 prim - human - lung-fibroblast
 prim - human - lymphocytes
 prim - human - PBMC
 prim - human - PBMC
 prim - human - peridontal fibroblast
 prim - human - skin, foreskin-fibroblast-1
 prim - human - skin, foreskin-fibroblast-2
 prim - mouse - bone marrow
 prim - mouse - embryonal stem cells-MPI-FR
 prim - mouse - fibroblast; perodontal ligament
 prim - mouse - MEF embryonal fibroblast
 prim - mouse - MEF embryonal fibroblast; transfected
 prim - mouse - spleen, MNC
 prim - mouse - spleen-isolated
 prim - mouse - spleen
 prim - mouse - thymocytes
 T cells - mouse - in lymphoma

TT - yeast

Yeast - Candida albicans
 Yeast - S cerevisiae Montpellier
 Yeast - Saccharomyces cerevisiae
 Yeast - Schizosaccharomyces pombe

CASY TTC | SPECIMEN MEASURED – SETUP FILES



TTC - Algae / Plankton / Protozoa

Alexandrium ostenfeldii - dino flagelate
 Blue algae - Cyano bacter synechocystis
 Cheatoceros calcitrans - fresh culture
 Cheatoceros calcitrans - old culture
 Chlorella vulgaris - fresh culture
 Chlorella vulgaris - old culture
 Desmodesmus
 Dictyostelium-viable
 Dino flagelates - fresh culture
 Dino flagelates - old culture
 Euglena gracilis - fixed in EtOH
 Euglena gracilis
 Gonyaulax polyedra
 Isochrysis galbana (Marine diatom)
 Ostreococcus tauri
 Phaeodactylum tricornutum
 Protozoa Dd - cultured in flask
 Protozoa Dd - fixed in EtOH
 Protozoa Tp - bioreactor - fast moving
 Protozoa Tp - bioreactor - immobilized
 Protozoa Tp - dead - fixed in EtOH
 Rhodomonas - mix viable-dead
 Rhodomonas - viable
 Scenedesmus subspicatus-01
 Scenedesmus subspicatus-02
 Tetracystis
 Tetrahymena thermophila - log phase
 Tetrahymena thermophila - stat phase

TTC - Bacteria

Acinetobacter baylyi DSM 14961
 Alcaligenes faecalis DSM 13975
 Aspergillus niger ATCC 6275
 Bacillus subtilis
 Bordetella petrii DSM 12804
 Chryso sporium merdarium DSM 62115
 Comamonas acidovorans
 Comamonas testosteronii DSM 1455
 E.coli - (1) - count
 E.coli - (2) - biovolume
 E.coli - K12 - wild type
 E.coli viable
 E.coli-BL21
 E.coli-SP6 - defect mutant
 Edwarsiella hoshinae DSM 13771
 Enterobacter cloacae DSM 16657
 Erwinia aphidicola DSM 19347
 Escherichia blattae DSM 4481
 Escherichia coli ATCC 11229
 Escherichia coli
 Klebsiella pneumoniae
 Kluyvera georgiana DSM 9409
 Lactobacillus - inaktivated
 Lactobacillus -5h in complete medium
 Lactobacillus
 Malikia spinosa
 Neisseria canis DSM 18000
 Pasteurella aerogenes DSM 10153
 Providenzia heimbachae DSM 3591

CASY TTC | SPECIMEN MEASURED – SETUP FILES

Prteus myxofaciens DSM 4482
 Pseudomonas aeruginosa
 Pseudomonas chloraphis subsp chloraphis DSM 50083
 Pseudomonas fluorescens DSM 6147
 Pseudomonas monteilii DSM 14164
 Pseudomonas mosselii DSM 17497
 Pseudomonas pseudoalcaligenes DSM 7521
 Pseudomonas putida
 Pseudomonas stutzerii
 Rahnella aquatilis DSM 4594
 Salmonella thyphimurium
 Serratia marcescens
 Staphylococcus aureus
 Staphylococcus capitis DSM 6180
 Staphylococcus carnosus DSM 11676
 Staphylococcus epidermidis
 Staphylococcus gallinarium DSM 206010
 Staphylococcus saprophyticus DSM 18669
 Staphylococcus simulans DSM 20322
 Staphylococcus xylosus DSM 6179
 Trichophyton mentagrophytes ATCC9533
 Yersina bercovieri DSM 18528

TTC - Blood Cells

Control Blood - human - after lysis
 Whole blood - dog-02
 Whole blood - dog-lysed-03
 Whole blood - dog-lysed
 Whole blood - dog
 B-lymphocytes - chicken
 Cord blood - human
 Granulocytes - basophilic – human
 Granulocytes - eosinophilic - human
 Macrophages - human
 Mono nuclear cells MNC - human
 Monocytes - human
 Neutrocytes - human
 Erythrocytes - chicken - early stage of differentiation
 Erythrocytes - human - membrane ghosts
 Whole blood, RBC – bovine
 Whole blood, RBC - human - high MCV
 Whole blood, RBC - human - low MCV

Whole blood, RBC - human - normal MCV
 Whole blood, RBC - pig
 Whole blood, RBC - sheep
 PLT - human - active
 PLT - human - mixture of inactive and active
 Platelets
 Red Blood Cells
 White Blood Cells

TTC - Blood parasites

Trypanosoma brucei
 Trypanosoma

TTC – Cell lines

293 (HEK) - human - embryonal kidney - mutant -E-cadherin
 293 (HEK) - human - embryonal kidney - wild type
 293 (HEK) - human - embryonal kidney-01-volume
 293 (HEK) - human - embryonal kidney-01
 293 (HEK) - human - embryonal kidney-02-volume
 293 (HEK) - human - embryonal kidney-02
 3T3 - mouse, swiss albino - fibroblast
 3T3-NIH - mouse, swiss albino - fibroblast - in EtOH
 3T3-NIH - mouse, swiss albino - fibroblast
 468 - mouse-mouse hybridoma - lymphoblast-01
 468 - mouse-mouse hybridoma - lymphoblast-02
 529 - mouse-mouse hybridoma - lymphoblast-01
 529 - mouse-mouse hybridoma - lymphoblast-02
 941 - mouse-mouse hybridoma - lymphoblast
 A-375 - human - skin, malignant melanoma
 A-431 - human - epidermoid carcinoma
 A549 - human - lung carcinoma - fixed in CASYblue
 A549 - human - lung carcinoma - mixture
 A549 - human - lung carcinoma – transfected
 A549 - human - lung carcinoma
 A9 - mouse - subcutaneous connective tissue-01
 A9 - mouse - subcutaneous connective tissue-02
 B9 - mouse - hybridoma - fixed in CASYblue
 B9 - mouse - hybridoma - fixed in EtOH
 B9 - mouse - hybridoma - viable-dead
 B9 - mouse - hybridoma – viable
 B9 - mouse - hybridoma 01
 B9 - mouse - hybridoma 02
 BHK-21 - hamster, syrian - kidney

CASY TTC | SPECIMEN MEASURED – SETUP FILES

BL-3 - bovine - lymphosarcoma - 01
 BL-3 - bovine - lymphosarcoma - 02
 BSC-1 - monkey, african green - kidney
 CCRF-CEM - human - T cell leukemia
 CHO - hamster, chinese - ovary cells - diploid
 CHO - hamster, chinese - ovary cells - serum free medium
 CHO - hamster, chinese - ovary cells - transfected
 CHO-K1 - hamster, chinese - ovary cells
 CHSE-214 - fish - Chinook salmon embryo cells-01
 CHSE-214 - fish - Chinook salmon embryo cells-02
 COS-7 - monkey, african green - kidney; fibroblast – dead
 COS-7 - monkey, african green - kidney; fibroblast - viable
 CTLL-2 - mouse - T-lymphocyte; cytotoxic - dead
 CTLL-2 - mouse - T-lymphocyte; cytotoxic - viable
 EOL-1 - human - acute myeloid leukemia
 FL - human - HeLa contaminant
 GOS-3 - human - astrocytoma; oligodendroglioma -fixed in CASYblue
 GOS-3 - human - astrocytoma; oligodendroglioma -viable-dead
 GOS-3 - human - astrocytoma; oligodendroglioma
 H4 - human - brain; neuroglioma - dead
 H4 - human - brain; neuroglioma - viable
 HeLa - human - cervix carcinoma
 Hep 3B - human - liver; hepatocellular carcinoma
 HEP-G2 - human - liver; hepatocellular carcinoma-fixed in CASYblue
 HEP-G2 - human - liver; hepatocellular carcinoma-viable-dead
 HEP-G2 - human - liver; hepatocellular carcinoma
 Hfob - human - foetal osteoblast
 HIGH 5 - insect - Trichoplusia ni -01
 HIGH 5 - insect - Trichoplusia ni -02 nuclei
 HL-60 - human - acute myeloid leukemia
 HT-29 - human - colon adenocarcinoma -fixed in CASYblue
 HT-29 - human - colon adenocarcinoma -viable-dead
 HT-29 - human - colon adenocarcinoma
 IMR-90 - human, caucasian - foetal lung fibroblast
 J82 - human - bladder carcinoma
 JURKAT - human - T cell leukemia-01
 JURKAT - human - T cell leukemia-02
 K-562 - human - chronic myeloid leukemia
 KC - insect - Drosophila embryo -01
 KC - insect - Drosophila embryo -02 nuclei
 L-1210 - mouse - lymphocytic leukemia
 L-929 - mouse - connective tissue fibroblast-01
 L-929 - mouse - connective tissue fibroblast-02
 L5178Y - mouse - lymphoma
 LNCaP - human - lymph node carcinoma
 M-NFS-60 - mouse - peripheral blood, virus induced myeloic leukemia
 MDA-MB486 - human - pleural effusion; adenocarcinoma - dead
 MDA-MB486 - human - pleural effusion; adenocarcinoma - viable-dead
 MDA-MB486 - human - pleural effusion; adenocarcinoma - viable
 MDBK - bovine - kidney cells
 MDCK - dog - kidney cells
 MEF-1 - mouse - embryo, fibroblast
 MGH - human - bladder carcinoma
 MOLT-4 - human - T-cell leukemia -fixed in EtOH
 MOLT-4 - human - T-cell leukemia -viable-dead
 MOLT-4 - human - T-cell leukemia
 MRC-5 - human - lung; fibroblast -fixed in EtOH
 MRC-5 - human - lung; fibroblast -mixture
 MRC-5 - human - lung; fibroblast
 MTC-M - mouse - thyroid carcinoma (1) aggr
 MTC-M - mouse - thyroid carcinoma (2) desintegration
 NAMALVA - human - Burkitt lymphoma
 P815 - mouse - mast cell; mastocytoma
 PER-C6 - human - embrional retinoblast - serum free medium
 PER-C6 - human - embrional retinoblast - treated with Accumax
 PER-C6 - human - embrional retinoblast
 RAJI - human - Burkitt lymphoma
 RAW 264.7 - mouse - monocyte - viable-dead
 RAW 264.7 - mouse - monocyte; macrophage; AML induced tumor
 RK13 - rabbit - kidney
 RPE-J - rat - eye (retina), pigmented epithelium
 RT-4 - human - urinary bladder transitional cell carcinoma

CASY TTC | SPECIMEN MEASURED – SETUP FILES

RTG-2 - fish, rainbow trout - gonadal cells
 SF-21 - insect - fall army worm -02 nuclei
 SF-21 - insect, fall army worm -01
 Sf9 - insect, fall army worm - infected, baculo virus
 Sf9 - insect, fall army worm - non infected
 Sf9 - insect, fall army worm -fixed in CASYblue
 SiHa - human - cervix; squamous cell carcinoma
 SP2-0 - mouse - myeloma - 01
 SP2-0 - mouse - myeloma - 02
 sp2-0-test-dead-01
 sp2-0-test-viable-01
 sp2-0-test-viable-dead-01
 TF-1 - human - erythroleukemia
 U-937 - human - histiocytic lymphoma -02
 U-937 - human - histiocytic lymphoma -camtothecin induced
 U-937 - human - histiocytic lymphoma -viable-dead
 U-937 - human - histiocytic lymphoma -viable
 U-937 - human - histiocytic lymphoma
 V-79 - hamster - lung fibroblast
 Vero - monkey, african green - kidney cells -fixed in EtOH
 Vero - monkey, african green - kidney cells
 VERO-B4 - monkey, african green - kidney cells
 Y79 - human - eye (retina), retinoblastoma -viable-dead
 Y79 - human - eye (retina), retinoblastoma
 YB2-0 - rat-rat - hybridoma -fixed in EtOH
 YB2-0 - rat-rat - hybridoma -viable-dead
 YB2-0 - rat-rat - hybridoma
 AG8 dead
 AG8 viable
 C168J - live cells
 C168J - fixed in CASYblue
 HUH7 - viable-dead
 HUH7 - viable
 LADA - viable
 MHH - dead-viable
 MHH - fixed in EtOH
 MHH - viable
 NS-1 - viable
 NSO - viable
 S180 - dead-viable
 S180 – dead
 S180 – viable

TTC - Fungi / Spores

Coniothyrium minitans
 Puccinia coronata

TTC - Primary Cells

Adipocytes - human
 Bone marrow - human - fixed in EtOH
 Bone marrow - human - viable cells
 Chondrocytes - human - 1. passage
 Chondrocytes - human - 3. passage
 Chondrocytes - human - fixed in EtOH
 Ehrlich ascites tumor cells - mouse
 Endothelial cells - mouse
 Endothelial cells 01 - human - fresh isolated
 Endothelial cells 02 - human - 1.passage
 Fibroblasts - human - inactiv
 Fibroblasts - human
 Glioma cells - mouse
 Hepatocytes - rat
 Hyalocytes - pig
 Keratinocytes - human - precursor cells
 Keratinocytes – human
 Keratinocytes - mouse
 Lung cells - human
 Lymphosarcoma cells - mouse
 Osteoblasts - human
 Pancreas cells - rat
 PBMC human EtOH fixiert
 PBMC human viable cells
 PBMC- human - isolated by Ficoll-01
 PBMC- human - isolated by Ficoll-02
 PBMC- human - lysed by saponin
 Renale proximale tubulus epithelial cells - human
 Spleen cells - mouse - fixed in CASYblue
 Spleen cells - mouse - freshly isolated
 Spleen cells, murine - mouse
 Thymocytes - mouse - freshly isolated
 Thymocytes - mouse -fixed in CASYblue
 Thymocytes - mouse -viable-dead
 Ziliar mussel cells - human

CASY TTC | SPECIMEN MEASURED – SETUP FILES

TTC - Sperms

sperms-Boar-01
sperms-Boar-02
sperms-Boar-03
sperms-Bull-01
sperms-Bull-02
sperms-Bull-03
sperms-Human-P01-hypoton
sperms-Human-P01-isoton
sperms-Human-P02-isoton
sperms-Human-P03-filtration
sperms-Human-P03-native
sperms-Human-P04-filtration
sperms-Human-P04-native
sperms-Stallion-01
sperms-Stallion-02

TTC - Stem Cells

Blood stem cells - human
ES DN 62
ES DN72 P29
stem cells - human, cord blood - dead and live
stem cells - human, cord blood - viable

TTC - Worm Eggs

worm eggs - trichuris suis-02
worm eggs - trichuris suis

TTC - Yeast

Candida albicans
Candida glabrata
Saccharomyces cerevisiae ATCC 1333
Saccharomyces cerevisiae DSM 70449
yeast - brewery - strain a
yeast - brewery - strain b
yeast - s. cerevisiae
yeast -01

CASY | MEET THE LEGEND



Questions? Remarks? Get in Touch!

We are looking very much forward to getting in contact with you.
Let's discuss your specific application / field of research!

- learn how CASY supports your daily laboratory cell-based work
- learn more about CASY models and CASY technology
- schedule a product demonstration
- request an offer

email: casy@ols-bio.de | fon: +49-421-276169-0

Here at OLS, we have experienced scientific staff happy to provide any assistance with your specific questions!

More resources

cellcounting.de is your online resource for your CASY. You'll find manuals, videos and a lot of support-topics. And even more: Tips, Tricks, new applications – and a direct contact for you to make your work an outright satisfying experience.

Watch CASY hands-on videos: bit.ly/cellcounting

Stay updated: register for our free newsletter: ols-bio.de/newsletter

www.cellcounting.de

Mail: info@ols-bio.de | www.ols-bio.de
fon +49-421 276169-0 | fax +49-421 276169-19
OLS OMNI Life Science GmbH & Co KG