


# CASY<sup>®</sup>-Technology

Applications Virus Infection of Cells

## Virus infection of Sf9 insect cells (*Spodoptera frugiperda*) infected with Baculo-Virus

With the standardised cell volume analysis of the CASY<sup>®</sup>-Technology it was for the first time possible to calculate the virus titer (MOI) by analysing the changes in cell volume of Sf9 cells infected with Baculo-Virus (insect cells – *Spodoptera frugiperda*).



By CASY<sup>®</sup> Model TTC using a 150 µm capillary

Copyright<sup>®</sup> by innovatis AG